28U11/28U23



28U22

28U11 WalkOn, 28U23 WalkOn Trimable


28U22 WalkOn Flex

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	13
IT Istruzioni per l'uso	18
ES Instrucciones de uso	22
PT Manual de utilização	27
SV Bruksanvisning	32
DA Brugsanvisning	36
NO Bruksanvisning	41
FI Käyttöohje	45
PL Instrukcja użytkownika	50
HU Használati utasítás	55
CS Návod k použití	59
RO Instrucțiuni de utilizare	64
HR Upute za uporabu	69
BG Инструкция за употреба	73
TR Kullanma talimatı	78
RU Руководство по применению	83
JA 取扱説明書	88
ZH 使用说明书	93

28U11* 28U22* 28U23*	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size	Höhe / Height	
	*=L36–39	links/left	36–39	35 cm	13.8 in
	*=R36–39	rechts/right	36–39	35 cm	13.8 in
	*=L39–42	links/left	39–42	36 cm	14.2 in
	*=R39–42	rechts/right	39–42	36 cm	14.2 in
	*=L42–45	links/left	42–45	37,5 cm	14.8 in
	*=R42–45	rechts/right	42–45	37,5 cm	14.8 in
	*=L45–48	links/left	45–48	39 cm	15.4 in
	*=R45–48	rechts/right	45–48	39 cm	15.4 in

Material	28U11 WalkOn, 28U23 WalkOn Trimable
	Carbonfaserverbundstoff/carbon fibre composite; Wadenband/calf strap: PCM® (Phase change material), ein temperaturregulierendes Material/a temperature regulative material, Mikroklett/micro hook and loop, weiches Klettband/soft velcro tape
	28U22 WalkOn Flex
	Glas-Carbonfaser-Gemisch/glass carbon fibre composite; Wadenband/calf strap: PCM® (Phase change material), ein temperaturregulierendes Material/a temperature regulative material, Mikroklett/micro hook and loop, weiches Klettband/soft velcro tape
	Latex-frei: Dieses Produkt enthält kein Naturgummi./ Latex-free: This product does not contain natural rubber.

INFORMATION

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INFORMATION

Last update: 2019-12-17

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.

- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 28U11 WalkOn, 28U23 WalkOn Trimable and 28U22 WalkOn Flex ankle-foot orthoses.

2 Intended use

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Drop foot, without or with mild spasticity

Caused for example by:

- Peroneal nerve injury
- Stroke
- Traumatic brain injury
- Multiple sclerosis
- Neuromuscular atrophy

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

- Moderate to severe spasticity in the lower leg
- Moderate to severe oedema
- Moderate to severe foot deformities
- Leg ulcers

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

2.4 Restrictions for use

The orthosis must not be used for the prevention of contractures (shortening) etc.

The orthosis must not be used for sports activities that involve jumping, sudden movements or rapid step sequences (e.g. basketball, badminton, riding sports). Sports activities should generally be discussed with the patient.

2.5 Effects

- Supports dorsiflexion
- Promotes dynamic rollover
- Promotes a physiological gait pattern
- Reduces the risk of falling

2.6 Additional options

- 28Z10 Lateral pronation strap

3 Safety

3.1 Explanation of warning symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

⚠ WARNING

Operation of motor vehicles

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

⚠ CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Use

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.

- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

4.2 Adaptation

Selecting footwear

NOTICE

Use without suitable footwear

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.
- ▶ Choose a sturdy lace-up shoe with solid heel cup to ensure the optimum effectiveness of the orthosis. The effective heel height should be between 0.5 cm and 1.5 cm.

Adapting the orthosis

⚠ CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

NOTICE

Damage due to improper handling (heat treatment, drilling holes)

Loss of orthosis functionality

- ▶ Do not carry out any heat treatment on the orthosis. The orthosis is not thermoformable.
 - ▶ Do not drill any holes in the bottom of the sole or the connecting elements as holes break the fibres and weaken the component.
- 1) Put the orthosis on the patient's foot and determine the relation of the longitudinal arch to the edge of the heel.
 - 2) If the insole in the patient's shoe can be removed, use it as a template to mark the correct size on the foot plate or transfer the contours of the foot directly to the plate with a pen (see fig. 1).
 - 3) If using the insole, pay attention to the position of the insole on the foot plate and the relation of the longitudinal arch to the edge of the heel. The orthosis should not come into contact with the medial malleolus, and the movement of the patient's Achilles tendon should not be restricted.
 - 4) **28U11:** Sand the foot plate up to the marking (see fig. 2).
28U22, 28U23: Use scissors to trim the foot plate along the marking (see fig. 2).
DANGER! When removing material to adjust the width of the foot plate, do not remove more than necessary on the lateral side to prevent the spring pressing against the medial malleolus. If the plate needs to be narrower, remove material on both sides.
 - 5) **28U11:** To adjust the length of the foot plate, focus on sanding the posterior section of the plate so that the spring is positioned behind the medial malleolus and pressure on the ankle and the medial malleolus is prevented (see fig. 3).
28U22, 28U23: To adjust the length, cut the foot plate on the posterior section of the plate (see fig. 3).
 - 6) **Optional:** If the lower leg shell needs to be adapted to the patient, sand the shell.

- 7) If the orthosis has sharp edges, smooth down the sanded edges.
- 8) **For foot deformities:** Correct the foot position with an insole or a custom moulded device.
INFORMATION: If the foot cannot be corrected by an insole, the orthosis and a sturdy shoe, the WalkOn should not be used.
- 9) Apply the textile component to the orthosis at calf height.
- 10) To adjust the length of the calf strap, put the orthosis on the patient.
- 11) **INFORMATION: The calf strap is the ideal length when it covers the hook-and-loop closure and has an excess length of 2 cm at most.**
Determine the length of the strap and shorten it with scissors, taking the patient's individual characteristics into account (see fig. 4).
- 12) Check that the orthosis fits correctly on the patient.

4.3 Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

- 1) Remove the insole from the shoe of the side to be fitted.
- 2) Open the calf strap.
- 3) Insert the orthosis into the shoe.
- 4) Put the orthosis on together with the shoe.
- 5) Fasten the calf strap.

4.4 Removal

- 1) Open the calf strap.
- 2) Remove the shoe together with the orthosis.
- 3) **Optional:** Remove the orthosis from the shoe and re-insert the insole of the shoe.

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

- ▶ Ottobock recommends cleaning the orthosis twice a week.

Textile component

- ▶ **NOTICE! Detergent residues may cause skin irritation and wear of material.**

- 1) Remove the textile component from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the textile component in warm water at approx. **40 °C** with a standard mild detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

Orthosis frame

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

7.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

7.3 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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