

50A10 Omo Immobil

50A11 Omo Immobil Rotation

DE	Gebrauchsanweisung	4
EN	Instructions for use	12
FR	Instructions d'utilisation	19
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	34
PT	Manual de utilização.....	41
NL	Gebruiksaanwijzing.....	49
SV	Bruksanvisning	56
DA	Brugsanvisning.....	63
FI	Käyttöohje.....	70
PL	Instrukcja użytkowania.....	77
HU	Használati utasítás	84
CS	Návod k použití.....	91
HR	Upute za uporabu	99
BG	Инструкция за употреба	105
TR	Kullanma talimatı.....	112
RU	Руководство по применению	119
JA	取扱説明書.....	127
ZH	使用说明书.....	134

	Größe / Size		Körpergröße / body height (m / inch)	
	50A10=S	50A11=L/R – S	~1.50 – 1.65	~58.5 – 64.5
	50A10=M	50A11=L/R – M	~1.65 – 1.75	~64.5 – 68.5
	50A10=L	50A11=L/R – L	From ~1.75	From ~ 68.5

L/R: left or right

Deutsch

INFORMATION

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INFORMATION

Last update: 2017-01-30

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

Explanation of symbols

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

1 Description

1.1 Indications for use

The 50A10 Omo Immobil/50A11 Omo Immobil Rotation shoulder joint brace is to be used **solely** for the orthotic fitting of the upper limb and is only suitable for contact with healthy skin.

1.2 Indications

1.2.1 50A10 Omo immobil

Post-operative and post-traumatic conditions of the shoulder and arm region that require abduction positioning of the arm, including conditions following rotator cuff surgery, reduction of posterior dislocation, subcapital fractures of humerus, shoulder impingement syndrome, shoulder replacement surgery or operative shoulder stabilization.

1.2.2 50A11 Omo immobil rotation

All indications that require the immobilization of the shoulder joint with an external rotation of 10°, 20° or 30°, including conditions following anterior dislocation or Bankart lesion (conservative treatment and post-operative immobilization). Indication must be determined by the physician.

1.3 Contraindications

1.3.1 Relative contraindications

In case of the following indications, consultation with a physician is required: Skin diseases/injuries, inflammatory symptoms, prominent scars on body parts to which the brace is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the brace is to be applied; sensory and circulatory disorders of the arm/hand.

1.4 Function

The 50A10 Omo Immobil shoulder joint brace allows for functional positioning of the shoulder through abduction by 0°, 15°, 30°, 60° or 90°. In this way, it effectively relieves the shoulder joint and directly supports the rehabilitation process. The 50A11 Omo Immobil Rotation also allows for immobilization with an external rotation of 10°, 20° or 30°. Any combination of abduction and external rotation angle can be used.

1.5 Scope of delivery and design (Fig. 1)

1.5.1 Scope of delivery for 50A10 Omo Immobil

- (1) Tube frame with fabric cover, waist belt, shoulder strap, shoulder pad and hand pad
- (2) Functional elements for setting abduction:
 - Double-sided hook strip: 0°
 - Small abduction wedge: 15°
 - Large abduction wedge: 30°, 60°, 90°
- (3) 2 straps with 2 double-sided hook strips for relieving pressure on the shoulder and neck region (optional)
- (4) Optional upper arm strap

1.5.2 Scope of delivery for 50A11 Omo Immobil rotation

- (1) – (4) identical to scope of delivery for 50A10 Omo Immobil
- (5) 29A222 Set for setting external rotation:
 - Element for 10° external rotation
 - Element for 20° external rotation
 - Element for 30° external rotation

1.5.3 Design

The 50A10 Omo Immobil/50A11 Omo Immobil Rotation shoulder joint brace is supplied ready for trial fitting on the right side.

The 50A10 Omo Immobil features a tube frame and a fabric-lined foam cover for padding purposes. Other functional elements, such as fixation straps, hand pad and abductor can be attached to the brace with hook and loop patches.

The 50A11 Omo Immobil Rotation brace also features various elements for setting the external rotation to 10°, 20° or 30°. When required, an element can be attached to the tube frame by undoing two locking pins.

2 General safety instructions

INFORMATION

The patient is to be instructed in the proper use/care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

CAUTION

Risk of injury as a result of improper use. The product is designed for use on **one** patient only. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

A brace/support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly.

NOTICE

Damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions.

3 Handling

3.1 Size selection

The size of the brace is selected according to body height.

3.2 Adapting and applying the product

The initial fitting and application of the product must be carried out by trained, qualified personnel only.

3.2.1 Adapting the product

Setting the body side

The shoulder joint brace comes already adjusted for the right arm. To set it up for the left side, remove the waist belt and the hand pad and attach them laterally reversed on the opposite side of the torso component (Fig. 2). When reattaching the waist belt, make sure that the strap guide loop is at the front (anterior). Remove the arm straps and undo the shoulder strap closure. Do not remove the shoulder pad from the shoulder strap. Once adjusted for the left shoulder, turn the pad on the strap so that seam is away from the body.

Setting the abduction angle

The desired abduction angle can be set with the abduction wedges that attach to the brace with hook and loop patches. To do this, unfold the brace and position the wedge at the upper pivot point of the joint (where the torso support and arm rest connect) on the arm rest (Fig. 3). Fold back the arm rest. If 0° abduction is desired, attach the double-sided hook strip at the position described above.

Shortening the shoulder strap or waist belt

INFORMATION To produce the desired effect, the brace must always be firmly applied. Custom adjustment is an important prerequisite for successful treatment.

If necessary, the waist belt and should strap can be shortened by an orthotist. To do this, remove the hook and loop closure from the strap, shorten the strap and re-attach the closure. The hook and loop closures should be located at the front (anterior side) of the torso for ease of use.

Setting external rotation on the 50A11 Omo Immobil Rotation

To mount the external rotation elements (**29A222**), remove the hook and loop closure that connects the upper and forearm pads so that the forearm pad can be pulled off the tube frame (Fig. 5). Remove the locking pins, then pull off the

forearm element and replace it with the desired external rotation element (Fig. 6). External rotation can be set to 10°, 20° or 30°, as desired. Fasten the tube frame using the locking pins and cover it with the forearm pad. Reconnect the forearm and upper arm pads with the hook and loop closure.

3.2.2 Applying the product

INFORMATION The brace should be applied and removed by two people in order to keep the patient's arm carefully and safely in place. The patient's elbow should be angled.

1. Open the shoulder strap closure as well as the closures of the waist belt and arm straps.
2. Position the brace under the arm, pull it up onto the torso and place the shoulder strap diagonally over the contralateral (healthy) side (Fig. 8). Position the shoulder pad (fastened to the shoulder strap) directly on the patient's neck (trapezius muscle). Guide the strap through the narrow return loop on the torso component and fasten the hook and loop closure (Fig. 8). The return loop is attached to the brace with a hook and loop patch that allows the position of the loop to be adjusted for individual adaptation to the patient.
3. Guide the waist belt through the return loop and close it (Fig. 9). The height of the waist belt that is attached to the brace with a hook and loop patch can be adjusted to the patient.
4. Correctly position the hand pad that is attached to the brace with a hook and loop patch as well (Fig. 9). Next, close the forearm straps, followed by the upper arm strap (Fig. 10). The direction for closing the arm straps can be chosen as desired (from the inside to the outside or from the outside to the inside). To ensure that the upper arm is held firmly in place in the brace, the included additional strap may also be used (Fig. 11).

INFORMATION If necessary (e.g., for adipose patients), the included additional straps can be used for relieving the pressure on the shoulder and neck region. Use the hook strips on the shoulder pad to fasten the strap. Cross the straps and attach them to the waist belt using the double-sided hook strips. Shorten the straps if necessary (Fig. 12).

Check that the brace is properly fitted before it is handed over to the patient for use.

INFORMATION Make sure that the brace does not cause pressure points in the underarm region and that the humerus head is not pushed upward.

If necessary, the orthotist can adjust the tube frame around the iliac crest for a more personal fit.

3.3 Instructions for use

NOTICE

Damage caused by wear and tear. Before each use, the brace must be checked for functional reliability and for possible wear or damage. If a closure or any other part of the brace shows any signs of wear, the brace should no longer be used.

Improper changes to the product are not permitted.

The brace must be used in strict conformity with the indication. The physician's instructions must be strictly observed.

Compared to a cast, the brace provides more comfort thanks to its light weight, hygienic features and adjustable arm position.

We recommend wearing a T-shirt or undershirt beneath the brace.

3.4 Care instructions

The fabric cover and straps/belt can be washed at 30 °C using a standard mild detergent. To do this, remove the hand pad and the abduction wedges from the brace. Then remove all hook and loop closures and well as the supporting tube frame. Soap residues can cause skin irritations. Therefore, rinse thoroughly after washing so that all residual detergent is removed. Do not expose to the sun and do not use oven or radiator heat for drying. The fabric cover can be dried in the dryer on a low setting (delicate cycle).

After drying, pull the fabric cover over the tube frame and reattach the hand pad, the hook and loop closures and the abduction wedges to the brace.

INFORMATION If necessary, a replacement cover can be ordered (article number 29A217). We recommend using the replacement cover while the fabric cover is being cleaned.

4 Further usage restrictions

CAUTION

Danger of injury due to improper use. Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the brace is used by another person.

The product is designed for use on **one** patient only. The daily duration of use and period of application are dependent on the medical indication.

Latex Safe: To our knowledge this product does not contain any natural rubber.

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

6.3 Local legal information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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Otto Bock HealthCare GmbH