
50R232 Smartspine LSO high


50R233 Smartspine LSO low

50R235 Smartspine LSO standard

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	12
IT Istruzioni per l'uso	16
ES Instrucciones de uso	20
PT Manual de utilização	24
NL Gebruiksaanwijzing	28
SV Bruksanvisning	32
DA Brugsanvisning	36
NO Bruksanvisning	40
FI Käyttöohje	44
IS Notkunarleiðbeiningar	48
PL Instrukcja użytkowania	52
HU Használati utasítás	56
CS Návod k použití	60
RO Instrucțiuni de utilizare	64
HR Upute za uporabu	68
SK Návod na používanie	72
BG Инструкция за употреба	76
TR Kullanma talimatı	80
JA 取扱説明書	84
ZH 使用说明书	88
KO 사용 설명서	91

	Größe / Size	Tailenumfang* / Waist circumference	
		cm	inch
	S	63,5–76,2	25–30
	M	76,2–88,9	30–35
	L	88,9–101,6	35–40
	XL	101,6–114,3	40–45
	XXL	114,3–127	45–50
	XXXL	127–139,7	50–55

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INFORMATION

Last update: 2018-05-08

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 50R232 Smartspine LSO high, 50R233 Smartspine LSO low and 50R235 Smartspine LSO standard lumbosacral braces.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

50R232 Smartspine LSO high

- Stenosis of the lumbar spine
- Degeneration of the lumbar intervertebral discs
- Spondylolisthesis of the lumbar spine
- Facet joint syndrome of the lumbar spine
- Spondylolysis of the lumbar spine
- Stable vertebral body fractures of the lumbar spine (with no neurological deficits)
- Postoperative stabilization after spinal fusion
- Low back/thoracic spine pain (mechanical)
- Compression fractures of the thoracic and lumbar spine (with no neurological deficits)
- Intervertebral disc prolapse in the thoracic and lumbar spine

50R233 Smartspine LSO low

- Stenosis of the lumbar spine
- Degeneration of the lumbar intervertebral discs
- Spondylolisthesis of the lumbar spine
- Facet joint syndrome of the lumbar spine

50R235 Smartspine LSO standard

- Stenosis of the lumbar spine
- Degeneration of the lumbar intervertebral discs
- Spondylolisthesis of the lumbar spine
- Facet joint syndrome of the lumbar spine
- Spondylolysis of the lumbar spine
- Stable vertebral body fractures of the lumbar spine (with no neurological deficits)

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the support is fitted, tumors, lymphatic congestion – including unclear soft tissue swelling distal to the support, and dysesthesia in the trunk.

2.4 Mechanism of Action

The brace promotes sensorimotor function and activates muscular stabilization of the lumbar spine. It supports pain relief and relieves the lumbar spine by increasing intra-abdominal pressure and by providing external stabilization.

The application of compression can be individually controlled with the patented "Mechanical Advantage Pulley System".

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE

Contact with ointments, lotions, or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to ointments, lotions, or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the waist circumference approx. **2 cm** above the iliac crest.
- 2) Determine the brace size (see size chart).

4.2 Adapt to patient

The height of the back plate can be customized to the needs of the patient.

- 1) Loosen the hook-and-loop closure of the back pad.
- 2) Mount the elastic bands of the abdominal strap at the desired position (see fig. 1).
- 3) Close the hook-and-loop closure of the back pad.

4.3 Application

⚠ CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the brace is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > The patient is standing.
- > The hook-and-loop closure of the abdominal strap is open.
- > The pull cable is fully loosened.
- 1) Center the back plate of the brace (see fig. 2).
- 2) Guide the brace around the body and apply equal force to both sides to pull it forward.
- 3) Attach the right side to the left side of the abdominal strap (see fig. 3). Ensure that the hook-and-loop sections lie flat against each other and are not in contact with clothing or skin.
- 4) Pull the end of the pull cable forward and attach it to the abdominal strap with the hook-and-loop fastener (see fig. 4).

- 5) **Optional – shortening the pull cable:** Pull the winding card out of the pocket (see fig. 5), wrap the pull cable around the winding card (see fig. 6), and put the winding card back into the pocket.

4.4 Optional – only 50R232 Smartspine LSO high: Remove components

The side parts, the abdomen part, and the back part can be removed. The removal of the components is directed by the physician.

Side parts

- 1) Open the hook-and-loop closures between the side parts and back part.
- 2) Remove the side parts. In order to do so, pull the side parts over the abdomen part (see fig. 7).

Abdomen part

- 1) Open the hook-and-loop closures on the abdomen part (see fig. 8).
- 2) Remove the abdomen part.

Back part

- 1) Open the hook-and-loop closure on the back part.
- 2) Remove the elastic bands from the back plate (see fig. 1).
- 3) Close the hook-and-loop closures on the back part.
- 4) Remove back part. In addition, pull the abdominal strap out of the back part.

4.5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Remove all boards/plates.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the brace in warm water at **30 °C** with standard mild detergent.
- 4) Rinse well.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

6.3 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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