ottobock.



Malleo TriStep 50S8

DE	Gebrauchsanweisung	6
EN	Instructions for Use	11
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	27
PT	Manual de utilização	33
NL	Gebruiksaanwijzing	38
SV	Bruksanvisning	44
DA	Brugsanvisning	49
NO	Bruksanvisning	54
FI	Käyttöohje	58
PL	Instrukcja użytkowania	63
HU	Használati utasítás	68
CS	Návod k použití	73
RO	Instrucțiuni de utilizare	79
TR	Kullanma talimatı	84
EL	Οδηγίες χρήσης	89
RU	Руководство по применению	95
JA	取扱説明書	101
ZH	使用说明书	106

Artikelnr. / Art. no.	Seite / Side	Größe / Size	Knöchelumfang cm / Circumference inch	
50S8=	L / R	- XS	25,5-29,0	10.0 - 11.5
50S8=	L / R	- S	29,5-31,5	11.6 - 12.5
50S8=	L / R	- M	32,0-34,0	12.6 - 13.5
50S8=	L / R	- L	34,5-36,5	13.6 - 14.5
50S8=	L / R	- XL	37,0-39,0	14.6 - 15.5

(A)

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Explanation of Symbols

English

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Date of the last update: 2013-03-18

- Please read this document carefully.
- Follow the safety instructions.

1 Description

INFORMATION

The instructions for use are to be handed over to the patient.

1.1 Intended Use

The 50S8 Malleo TriStep Ankle Orthosis is to be used **exclusively** for the orthotic fitting of the ankle joint and is only suitable for contact with healthy skin.

1.2 Indications and Effects

- Acute and chronic capsular ligament instabilities, ligament insufficiencies and injuries of the ankle joints (upper ankle joint / lower ankle joint) as well as severe distortions and lateral ligament ruptures
- Post-operative protection after seam or reconstruction of the lateral ligaments of the upper ankle joint
- Diseases and injuries of the peroneal tendons
- Inflammatory and degenerative joint diseases of the ankle joints (upper ankle joint / lower ankle joint) with recurring irritations or chronic reactive effusions

Specific indications must be determined by the physician.

The fitting takes place in 3 stages:

Stage I – Inflammatory stage (Fig. 1): The foot is immobilized in the neutral position (90°) with slight pronation using the foot shell and the pronation strap. Additional stability is achieved through the autoadaptive stabilizing elements on the sides. In this stage, the 50S8 Malleo TriStep cannot be worn within a shoe.

Stage II – Proliferation stage (Fig. 2): After the initial function as an immobilization orthosis in stage I, the foot shell can be removed during the day. With controlled weight bearing, mobilization / activation comes to the fore in stage II. For stabilization, the 50S8 Malleo TriStep is worn with stabilizing strap and with the autoadaptive stabilizing elements on the sides. At night, the foot shell is adapted to immobilize the ankle joint and to prevent recurrence.

Stage III – Remodulation stage (Fig. 3): At the end of the healing process, the stabilizing elements on the sides and the stabilizing strap are removed, and the 50S8 Malleo TriStep can be used as support and to prevent new injuries, e.g. during sports or other activities. As an option, the stabilizing elements on the sides and the stabilizing strap can be attached, as necessary.

1.3 Contraindications

1.3.1 Absolute Contraindication

The orthosis is not suitable for use as a dorsiflexion-assist orthosis in case of peroneal paralysis.

1.3.2 Relative Contraindications

In case of the following indications, consultation with a physician is required: Skin diseases /injuries, inflammatory symptoms, prominent scars on body parts to which the orthosis is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the orthosis is to be applied; sensory and circulatory disorders of the leg.

1.4 General Safety Instructions

INFORMATION

The patient is to be instructed in the proper use /care of the product. The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

⚠ CAUTION

Risk of injury as a result of improper use. The product is designed for use on **one** patient only. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

An orthosis / support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep it away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions.

1.5 Scope of Delivery and Design

The 50S8 Malleo TriStep Ankle Orthosis is supplied ready for stage I use. The velour fabric material consists of skin-friendly polyamide. The integrated plastic insert (polyethylene) of the orthosis is anatomically pre-shaped. If necessary, the orthotist may slightly modify its shape by cold forming. The

individually adjustable hook and loop closures are arranged to close in opposite directions to prevent the orthosis from slipping.

2 Handling

2.1 Size Selection (Fig. A)

Size selection of the orthosis is based on the ankle joint circumference.

2.2 Adapting and Applying the Product

The initial fitting and application of the 50S8 Malleo TriStep Ankle Orthosis must be carried out by qualified personnel. On delivery, the Malleo TriStep orthosis must be tested for the correct fit.

Stage I – Inflammatory stage: To apply the orthosis first open the cross-wise placed straps and then the individual hook and loop closures completely. This allows to comfortably get into the orthosis. Place your foot in the orthosis (Fig. 4) and first close the middle hook and loop closure by threading the strap through the loop and closing it accordingly (Fig. 5). Then close the two other hook and loop closures (Fig. 6). Continue with applying at first the outer (lateral) cross strap (Fig. 7) and then the inner (medial) cross strap (Fig. 8). Finally apply the pronation strap starting with attaching the lower part to the orthosis enabling you to pull the foot with the upper part of the strap into pronation (Fig. 9).

Stage II – Proliferation stage: Mobilization / activation comes to the fore in this stage. The orthosis is worn with the stabilizing strap and stabilizing elements. First, open the pronation strap and the cross straps and remove the foot positioning splint from the orthosis (Fig. 10). Then attach the stabilizing strap to the orthosis from the rear (Fig. 11). The "Ottobock" print must be readable (Fig. 11). Moreover, the red dots on the lateral pockets should precisely fit the holes in the stabilizing strap. To close the stabilizing strap place the outer (lateral) strap over the instep, lead it over the inner edge of the foot and then along beneath the foot (Fig. 12) and attach it to the lateral ankle (Fig. 13). Proceed laterally reversed with the strap on the other side (Fig. 14).

In the course of the healing phase the more stable, red stabilizing elements on the sides can be replaced by the more flexible, grey elements.

The point in time when this exchange can take place has to be determined by the treating physician.

Stage III – Remodulation stage: The stabilizing elements can be removed completely (Fig. 15), and the 50S8 Malleo TriStep is used as a soft orthosis in this therapy stage. Depending on the patient's feeling of stability and corresponding recommendation by the physician, the stabilizing strap and stabilizing elements on the sides can additionally be applied.

2.3 Instructions for Use

NOTICE

Damage caused by wear and tear. Before each use, the orthosis must be checked for functional reliability and for possible wear or damage. If a closure or any other part of the orthosis show any signs of wear, the orthosis should no longer be used. Improper changes to the product are not permitted.

The orthosis must be used in strict conformity with the indication. The physician's instructions must be strictly observed.

2.4 Care Instructions (Fig. B)

Hand-wash in 30°C warm water with a standard mild detergent. Close the hook and loop closures prior to washing. Soap residues can cause skin irritations. Therefore, rinse thoroughly so that all residual detergent is removed. Let air dry. Do not expose to the sun and do not use oven or radiator heat for drying.

The foot shell can be cleaned with a damp piece of cloth.

NOTICE

Damage caused by improper cleaning agents. The patient must be instructed not to clean the foot shell with any aggressive media.

3 Further Usage Restrictions

The product is designed for use on **one** patient only. The daily duration of use and period of application are dependent on the medical indication.

⚠ CAUTION

Danger due to improper use. Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the orthosis is used by another person.

FDA: This product / the package of the product contains synthetic latex.

4 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

5 CE conformity

This product meets the requirements of the 93 / 42 / EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH